The logo consists of a stylized, light blue graphic on the left, resembling a map or a vacuum seal, followed by the text "BodyMap" in a light blue, sans-serif font.

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Manufacturer takes responsibility for proper usage of the product only when the product was purchased from qualified REH4MAT employee or at specialized medical store.

REH4MAT reserve the right to introduce technical and trading changes in the content of the instruction without warning.

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Product marked with the CE sign is Class I medical device, of rule 1, consistent with the requirements of Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices.

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REH4MAT Company introduced quality management system ISO 13485 regarding production, sales and service. The system is certified by Det Norske Veritas.

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
**BodyMap<sup>®</sup>** System refers to first class medical products which are designed for adults as well as children with cerebral palsy. All these products provide perfect stabilization during seating or lying position.

What is more, it is a perfect choice for the disabled due to the fact that it prevents bed sores.

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**BodyMap<sup>®</sup>** System uses vacuum method. This method refers to pillows and mattresses which are filled with polystyrene granules thanks to which products can model according to patient's body. Polystyrene granules move freely inside the device thanks to which it is possible to model the cushion. The cushion perfectly models according to patient's body when patient either sits or leans on the cushion. The shape of cushion doesn't change easily because of the fact that the air is sucked out of the cushion.

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
### VACUM VALVE

It is used to pump the air out of the cushion. Equipped with a chain that prevents the valve from being lost.

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### MAINTENANCE VALVE

Together with the granule bag makes it possible to control the amount of the granules in the BodyMap<sup>®</sup> cushions. The valve is a standard equipment in the selected BodyMap<sup>®</sup> cushions. In other cushions, the valve is made on request.

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### SELF-ADHESIVE FUNCTION

The bottom of selected cushions is covered with self-adhesive function. Various elements can be easily attached to it with Velcro tapes. In other cushions, the self-adhesive function is made on request.

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### TARGET GROUP

People requiring stabilization and compensation of deficits resulting from deformation of the osteoarticular system and weakening of postural muscles.

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### INDICATIONS FOR USE

- Neurological and orthopedic diseases requiring stabilization of the body or its parts.
- Conditions after surgical or orthopedic procedures requiring stabilization of the body or its parts.
- Posture defects.
- Muscle weakness Cerebral Palsy - various forms.
- Postural therapy.

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### CONTRAINDICATIONS FOR USE

- Allergic to neoprene.
- Incorrect positioning of a person using the BodyMap<sup>®</sup> System may lead to deformation of the joints, and keeping the person lying in one position for a long time significantly increases the risk of pressure ulcers.

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
Concerning shaping the cushion, it is highly important to open valve so as to let the air in. This activity is done in order to equalize pressure inside the device with atmospheric pressure.

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
In order to open valve it is necessary to put thumb on frontal valve surface. Moreover, use your middle as well as index finger to press blue ring (1) and simultaneously slip off valve from pipe (2).

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Be careful with the valve because it is small and it can be easily lost. Be cautious especially when you perform all activities outside.

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
Cushion without valve should be left for few minutes only so as to equalize pressures.

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Next, shape the cushion so that it fits to the user's body. In order to shape the cushion, spread the granules evenly and let the air out with a pump so that the cushion is still soft and it is possible to shape it. Next you can shape it by sitting or leaning on it - it depends on the type of the cushion. Now you can start to properly pump the air out of the cushion.

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
It is not allowed to press and load the BodyMap<sup>®</sup> device in case the air is inside and the vacuum valve is assembled. This may damage the pillow!

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Move piston of the pump forward and backward so as to let out the air.

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
Air should be let out until it is properly hard. The shape of cushion can be easily shaped during this activity. Pump out longer so as to make cushion harder and more solid.

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Once the shape as well as hardness of cushion is done, it is necessary to disconnect vacuum pump. Valve will automatically close air supply for cushion.

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### BM\_001 SMALL HAND VACUUM PUMP

Handy, small vacuum pump. Perfect for small and medium cushions. Each cushion includes this vacuum pump free of charge.

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### BM\_002 LARGE HAND PUMP

Very efficient manual vacuum pump. Perfect for mattresses and big cushions.

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### BM\_003 ELECTRIC PUMP

Handy and very efficient electric vacuum pump. Equipped with electric charger. Can work with AA battery also (not included).

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### BM\_005 ELECTRIC DUAL-FUNCTION PUMP

Electric dual-function pump allows precise pumping out or in the air in order to achieve the desired hardness. Strongly recommended for large cushions and mattresses.

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### BM\_006 AIR VACUUM BLOW GUN

The Air Vacuum Blow Gun helps to quickly increase the amount of granules in the cushion. It's especially recommended for large cushions and mattresses.

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**ATTENTION!**

To avoid falling down the granules during molding the cushions with backrest (B, C, AB, AC, etc.) please follow the steps below:

- place the cushion on flat surface,
- spread the granules inside the cushion evenly,
- vacuum the cushion a little bit so it can still be floppy and moldable,
- put the cushion on its place,
- start the final molding seating on it or mold it manually.

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
**ATTENTION!**

Applies to BodyMap<sup>®</sup> K pillow: Before cleaning the cover with VISmemo<sup>™</sup> thermal foam, remove the foam from the inner pocket of the cover. The foam can be wiped with a dampened, clean cloth. Subjecting it to improper cleaning may cause it to lose its properties or become damaged.

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**Other BodyMap<sup>®</sup> pillows:** The cover with VISmemo<sup>™</sup> thermoactive foam should be pulled off the pillow and washed by hand at up to 40°C, using a mild detergent, kneading it gently. Do not use fabric softener. Remove excess water by gently kneading the water (folding the cover may be helpful).

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**ATTENTION!**

VISmemo™ thermal foam is very delicate, so special precautions are necessary. Products with VISmemo™ thermoactive foam must not be washed in a washing machine, spun, wrung out.

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**ATTENTION!**

Regarding a waterproof cover with the thermo-active foam, it is not possible to remove the foam. In this case, the cover can be disinfected with the means shown in the table below.

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After pumping the air out, the cushion should keep its shape for at least 48 hours\*. In case it loses its shape sooner, check the cushion yourself by carrying out **THE AIR LEAKAGE TEST**:

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### **STEP 1. Check if the cushion has no visible mechanical damage, tearing or signs of wear and if the connective elements have not been damaged.**

If the cushion has a small tearing, you can try to fix it yourself using the REPAIR KIT, which is added to the cushion (follow the instructions included). If the damage is significant, contact the manufacturer to have the cushion repaired (warranty does not cover mechanical damage or tearing of the cushion).

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### **STEP 2. Carry out the maintenance of the tube, which is used for pumping the air out:**

- take the blue tube out of the cushion by pressing the blue clenching elements and remove the valve,
- cut off about 10 mm of the tube on both sides,
- insert the tube back into the cushion and put the valve on it.

After STEP 2, check again if the cushion is airtight by shaping the cushion and pumping the air out of it. If the cushion still does not keep the shape, proceed to STEP 3.

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\*There are a lot of factors which affect the shape of the BodyMap<sup>®</sup> cushion, such as the kind of the cushion, its size, frequency of use, ambient temperature and the weight of the patient. The cushion should keep its shape from several days to several weeks. If the cushion is not used and does not keep its shape in the period of 48 hours, it might mean that it is damaged. To check that, carry out the AIR LEAKAGE TEST which is described above. If the test confirms that the cushion is damaged, inform the manufacturer about that fact in order to have the cushion repaired or replaced.

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**STEP 3. If the cushion is equipped with the service valve, check if the valve is airtight:**

- open the valve cap with the use of a special key,
- check if the surface around the valve or the seal are not damaged,
- clean and wash the seal with soap or dishwashing liquid; put on and close the valve cap, paying particular attention that the granules do not get under the seal.

All of the abovementioned steps should be done carefully because the sudden or rapid movements may result in the granules pouring out of the cushion. Do not press the cushion when the valve cap is removed.

After STEP 3, check again if the cushion is airtight by shaping the cushion and pumping the air out of it. If the cushion still does not keep the shape, send the cushion back to the manufacturer to have it repaired.

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Maintenance valve is an innovative as well as easy solution used in BodyMap® system. It makes it possible to control the amount of granules in BodyMap® cushions. It is highly advantageous because then it is possible to adjust cushion according to your needs.

Valve is situated at the back of cushion thanks to which it is not very visible and, moreover, it is not an obstacle while using the product. What is more, it is really easy to use it – you just need to unscrew cork and either pour in or take out specific amount of granules. Additional amount of granules is packed in airtight bags that are available in two sizes. Remember that the size of bag depends on the size of BodyMap® cushion.

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Please follow these steps in order to pour in / take out granules:

1. Vacuum the cushion using the vacuum pump.

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2. Unscrew carefully valve's cork.

**CAUTION!** When valve is unscrewed it is forbidden to move pillow violently or press it because then granules may spill out.

There is a key added to the bag which facilitates opening and closing the service valve in order to add or decrease the amount of the granules

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3. Take out cork which secures funnel from bag with granules. Be careful so as not to spill out its contents.

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4. Insert tightly bag's funnel into hole from cushion at a depth of 1cm. Pull bag upwards, pour granules in as much as it is needed. If you want to take out some granules from cushion, it is necessary to insert funnel into hole from cushion and then turn cushion in a way that valve is placed down. Next, take out granules as much as it is needed.

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Be careful when screwing the service plug so that the granulate does not get under the gasket - it may cause the pillow does not keep the shape.

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Over time, the polystyrene granules, which fill the BodyMap cushions, may give signs of flattening. As a result, the volume of the cushions will decrease, which is their natural feature. If this situation arises, we recommend replacing or refilling the granules in the cushion. It can be easily done with a use of a service valve which most BodyMap<sup>®</sup> cushions are equipped with. In particular, all seats, mattresses and beans.

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In the cushions not equipped with a service valve it is not possible to refill the granules. We recommend you to consider the option of adding such a valve while ordering. The decreasing volume of the cushions is their natural feature and it shall not give rise to a complaint.

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Detailed information about the cushions equipped with a service valve can be found in the catalogue and on the website [www.bodymapsystem.com](http://www.bodymapsystem.com).

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Before using the BodyMap<sup>®</sup> products, always make sure that the place of use (e.g. bed, mattress, seat) is free from sharp edges, metal parts, springs and other elements which may damage the product. The manufacturer is not liable for the mechanical damage of the product!

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First, place the BodyMap<sup>®</sup> cushion on the chair, stroller etc.

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In order to minimize the risk of bedsores, we recommend using a cotton cover along with each BodyMap<sup>®</sup> cushion if the permanent duration of use is longer than 2 hours. For people using BodyMap<sup>®</sup> cushions permanently for more than 8 hours it is recommended to use it along with thermoelastic ViSmemo<sup>™</sup> cover. Caution should be taken in subjects to people with known bedsores - prone. In such a case it is strictly recommended to you use BodyMap<sup>®</sup> cushion with aeration-massaging pad for people who use BodyMap<sup>®</sup> cushions non-stop longer than 8 hours! Under any circumstances BodyMap<sup>®</sup> product cannot be fold, when the air is pump out and the vacuum valve is assembled.

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Product works correctly only when the device is adjusted according to patient's height as well as weight. Correct adjustment is the most important due to the fact that only then the product is safe. It is forbidden to force child to keep specific position because he/she will feel uncomfortable and it may lead to wrong posture.

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If the symptoms occur (for example skin allergies, intertrigo, bedsores), stop using the product immediately, take the user out of the product and lay them down. If the undesirable symptoms persist, contact a physician immediately. If the symptoms have disappeared, inform the attending physician about this fact in order to make a decision about the further rehabilitation process.

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Serious incidents involving this device should be reported both to manufacturer and to the competent authority in your country.

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It is not allowed to press and load the BodyMap<sup>®</sup> device in case the air is inside and the vacuum valve is assembled. In order to shape the cushion, there is necessary to spread the granulate evenly and to pump the air out slightly. By maintaining the cushion soft, you can shape the device according to your needs.

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
The Velcro surface of BodyMap cushions is used for positioning of cushions on the product or fix the supporting elements on it.

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There are a lot of factors which affect the shape of the BodyMap<sup>®</sup> cushion, such as the kind of the cushion, its size, frequency of use, ambient temperature and the weight of the patient. The cushion should keep its shape from several days to several weeks. If the cushion is not used and does not keep its shape in the period of 48 hours, it might mean that it is damaged and the manufacturer should be informed about that fact in order to repair or replace the faulty cushion.

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
## USER MANUAL

VACUUM POSITIONING SYSTEM BODYMAP FOR TRUNK, HEAD AND LIMBS STABILIZATION.

28/U/14

WERSJA / VERSION: 07.06.2023